

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, ORANGE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, ORANGE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043714
- Mfr. No.: E10-3021B-ORG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670178

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Rigid TQ Case for CAT with TekLok Safety Instruction Guide](#)
- [About Us](#)

Rigid TQ Case for CAT with TekLok Safety Instruction Guide

Introduction

Thank you for choosing the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. This guide provides essential safety instructions to ensure the safe and effective use of your product. Please read this guide carefully before using the case to ensure optimal performance and safety.

General Safety Guidelines

- Ensure the Rigid TQ Case® is used only for its intended purpose: to securely hold and provide rapid access to the CAT® Tourniquet.
- Regularly inspect the case for any signs of damage or wear. A damaged case may not provide adequate protection or functionality.
- Keep the tourniquet and case clean and dry to ensure proper operation.
- Store the case in a cool, dry place away from direct sunlight to maintain the integrity of the materials.
- For safety, always follow the manufacturer's guidelines for the use of the CAT® Tourniquet.

Specific Safety Precautions for Use

- Do not overload the case with items other than the CAT® Tourniquet as this may compromise its function.
- Ensure that the case is properly secured to your belt or vest to prevent accidental detachment during use.
- Avoid using the case in extreme environmental conditions that may affect its performance, such as extreme heat or cold.
- If you are using the case in a highstress situation, ensure you have practiced deploying the tourniquet to familiarize yourself with the operation.
- Keep the case and tourniquet out of reach of children to prevent accidental use.

Instructions for Installation and Usage

1. Mounting the Case:

- Choose a suitable location on your duty belt or MOLLE/PALS compatible vest.
- Use the appropriate mounting hardware to secure the case in place. Ensure it is tight and secure.
- For vertical, horizontal, canted, or inverted mounting, adjust the position according to your preference and comfort.

2. Accessing the Tourniquet:

- To deploy the tourniquet, grasp the exposed portion of the tourniquet firmly.
- With a swift pull, release the tourniquet from the case. The design allows for onehanded deployment even under challenging conditions.
- Practice this action regularly to ensure efficiency in an emergency situation.

3. Maintenance:

- Clean the case regularly with a damp cloth to remove dirt and debris.
- Check the mounting hardware periodically to ensure it remains secure.
- Replace any worn or damaged components immediately to maintain safety.

Disposal Instructions

- When the Rigid TQ Case® reaches the end of its life or is no longer needed, dispose of it in accordance with local regulations.
- Do not dispose of the case in regular household waste if it contains any hazardous materials.
- Consider recycling the materials if possible, following your local recycling guidelines.

Contact Information for Further Support

For questions regarding safety, usage, or maintenance, please refer to the manufacturer's website or contact customer support for assistance.

By following these safety instructions, you can ensure the safe and effective use of your Eleven 10® Rigid TQ Case® for the CAT® Gen 7. Always prioritize safety and training in the use of your equipment.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk