

# RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



## Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, MULTICAM](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043715
- Mfr. No.: E103021BMTC
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615672028

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Rigid TQ Case for CAT with TEKLOK Safety Instructions](#)
- [About Us](#)

# Rigid TQ Case for CAT with TEKLOK Safety Instructions

## Introduction

Thank you for choosing the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. This product is designed to provide easy access to your tourniquet in critical situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

## General Safety Guidelines

- Ensure that all users are familiar with the operation of the Rigid TQ Case®.
- Always inspect the case before use for any signs of damage or wear.
- Keep the case clean and free from debris to ensure proper function.
- Store the case in a cool, dry place when not in use.
- Do not modify the case in any way that may compromise its integrity or functionality.
- Follow all local laws and regulations regarding the use of tourniquets and first aid supplies.

## Specific Safety Precautions for Use

- The Rigid TQ Case® is designed for use by trained personnel. Ensure that you have received proper training in the application of tourniquets.
- Avoid exposing the case to extreme temperatures, chemicals, or direct sunlight for prolonged periods.
- Do not use the case if it shows any signs of damage, such as cracks or broken components.
- Keep the case out of reach of children and vulnerable individuals to prevent accidental use.
- Ensure that the tourniquet is correctly positioned within the case before securing it.

## Instructions for Installation and Usage

### Installation

1. Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
2. Use the provided mounting hardware to secure the Rigid TQ Case® in the desired position.
3. Ensure that the case is firmly attached and does not move during use.
4. If applicable, utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.

### Usage

1. To access the tourniquet, locate the exposed portion of the tourniquet in the case.
2. With one swift pull, remove the tourniquet from the case.
3. Follow established protocols for applying the tourniquet as needed.
4. After use, clean the case and inspect it for any damage before storing it again.

## Disposal Instructions

- Dispose of the Rigid TQ Case® in accordance with local regulations regarding plastic and medical waste.
- Do not incinerate the case, as it may release harmful fumes.
- If the case is damaged beyond repair, consider recycling the materials if possible.

## **Contact Information for Further Support**

For any inquiries or support regarding the Rigid TQ Case®, please refer to the manufacturer's contact details provided with the product. Ensure you have the product information ready for the most efficient assistance.

By following these guidelines, you can ensure the safe and effective use of the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. Always prioritize safety and proper training in first aid applications.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)