

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, BLACK](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043716
- Mfr. No.: E10-3021M-BLK
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670192

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Rigid TQ Case for CAT with TekLok Safety Instruction Guide](#)
- [About Us](#)

Rigid TQ Case for CAT with TekLok Safety Instruction Guide

Introduction

Thank you for choosing the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. This product is designed to provide secure and rapid access to your tourniquet in critical situations. Please read this safety instruction guide carefully to ensure safe use and proper handling of the product.

General Safety Guidelines

- Ensure that the Rigid TQ Case is used only for its intended purpose—holding and deploying a tourniquet.
- Regularly inspect the case for wear and tear. If damaged, discontinue use immediately.
- Keep the product out of reach of children to prevent accidental use or injury.
- Always follow proper first aid procedures when using a tourniquet.
- Report any unsafe product conditions or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Avoid exposing the Rigid TQ Case to extreme temperatures that may affect its integrity.
- Ensure that the tourniquet is properly secured within the case before use.
- Do not attempt to modify the Rigid TQ Case or its components.
- Use the product in accordance with the manufacturer's instructions to ensure safe operation.
- When deploying the tourniquet, ensure that your hands are clean and dry to maintain grip and control.

Instructions for Installation and Usage

Installation

1. Mounting Options:

- The Rigid TQ Case can be mounted on a duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- Choose a mounting position that allows easy access to the tourniquet.

2. Secure Attachment:

- Use the included Vibratite to secure all hardware in place once the desired configuration is obtained.
- Ensure that the case is firmly attached and does not move during use.

Usage

1. Accessing the Tourniquet:

- To deploy the tourniquet, grasp the exposed portion of the tourniquet.
- Pull swiftly to release the tourniquet from the case.

2. Deployment:

- Follow standard first aid procedures for applying a tourniquet.
- Ensure that the tourniquet is applied correctly and securely to control bleeding.

Disposal Instructions

- Dispose of the Rigid TQ Case in accordance with local regulations for plastic and electronic waste.
- Do not incinerate the product as it may release harmful substances.
- If the product is no longer usable, ensure that it is disposed of safely to prevent accidental use.

Contact Information for Further Support

For any inquiries or concerns regarding the Rigid TQ Case, please contact the manufacturer directly. Ensure to have your product details ready for assistance.

By following these guidelines, you can ensure the safe and effective use of your Eleven 10® Rigid TQ Case® for the CAT® Gen 7. Thank you for prioritizing safety and preparedness.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk