

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, COYOTE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043717
- Mfr. No.: E10-3021M-CYT
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670215

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [RIGID TQ Case for CAT with TEKLOK Safety Instructions](#)
- [About Us](#)

RIGID TQ Case for CAT with TEKLOK Safety Instructions

Introduction

Thank you for choosing the Eleven 10® RIGID TQ Case® for the CAT® Gen 7. This hardsided tourniquet case is designed for quick access to your tourniquet in critical situations. It is essential to use this product safely and correctly to ensure its effectiveness. This safety instruction guide provides you with important information regarding the proper use, installation, and disposal of the RIGID TQ Case®.

General Safety Guidelines

- Ensure that the product is used only for its intended purpose.
- Keep the RIGID TQ Case® out of reach of children and vulnerable individuals.
- Regularly check the condition of the case and the tourniquet to ensure they are in good working order.
- Report any unsafe conditions or accidents to the appropriate authorities.
- Stay informed about product recalls or safety updates through the EU Safety Gate platform.

Specific Safety Precautions for Use

- Verify that the tourniquet is compatible with the RIGID TQ Case® (CAT® Gen 7 and older generations).
- Ensure that the case is securely mounted to your duty belt, vest, or holster before use.
- Do not overload the case with additional items; it is designed solely for the tourniquet.
- Use the tourniquet only as instructed in the manufacturer's guidelines.
- In case of blood or wet conditions, practice using the case to ensure familiarity with onehanded deployment.

Instructions for Installation and Usage

Installation

1. **Select Mounting Position:** Choose a suitable location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
2. **Mounting the Case:**
 - For Police Duty Belt: Attach the BladeTech TEKLOK Belt Attachment to belts up to 2.25 inches.
 - For Safariland® Drop Leg Shroud: Position the case in the accessory position (mounting hardware sold separately).
3. **Secure the Hardware:** Utilize the included Vibratite to secure the mounting hardware in place once the desired configuration is obtained. Follow the instructions provided for proper installation.

Usage

1. **Accessing the Tourniquet:**
 - To deploy the tourniquet, grasp the exposed portion of the tourniquet firmly.
 - With a swift pull, the tourniquet will release from the case.
2. **Deployment:** Practice deploying the tourniquet in a safe environment to ensure proficiency, especially under stressful conditions.

Disposal Instructions

- Dispose of the RIGID TQ Case® and any associated components in accordance with local regulations.
- If the product is no longer usable, consider recycling materials where possible.
- Ensure that any used tourniquets are disposed of safely to prevent accidental use.

Contact Information for Further Support

For any questions or concerns regarding the RIGID TQ Case® or its usage, please refer to the manufacturer's website or customer service for assistance. Always ensure that you have the latest information regarding the product and its safety features.

Thank you for prioritizing safety and for your commitment to responsible usage of the RIGID TQ Case®.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk