

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, ORANGE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, ORANGE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043719
- Mfr. No.: E10-3021M-ORG
- Size: 1 Person
- Delivery weight: 0.088kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670239

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [RIGID TQ Case for CAT with TEKLOK: Safety Instruction Guide](#)
- [About Us](#)

RIGID TQ Case for CAT with TEKLOK: Safety Instruction Guide

Introduction

Thank you for choosing the Eleven 10 RIGID TQ Case for the CAT Gen 7. This hardsided tourniquet case is designed for rapid access to your tourniquet in critical situations. Please read this safety instruction guide carefully to ensure the safe and effective use of the product.

General Safety Guidelines

- Ensure that the RIGID TQ Case is used for its intended purpose only.
- Always inspect the case and the tourniquet before use to ensure they are in good condition.
- Keep the case and tourniquet out of reach of children when not in use.
- Follow local regulations and guidelines regarding the use of tourniquets and first aid supplies.
- Report any unsafe products or accidents to the appropriate authorities.

Specific Safety Precautions for Use

- **Proper Handling:** Always handle the tourniquet carefully to avoid injury. The case is designed for easy access, but ensure that your hands are dry and clean when deploying the tourniquet.
- **Mounting Safety:** When mounting the case, ensure that it is securely attached to your duty belt or vest to prevent accidental detachment during use.
- **Environmental Considerations:** Avoid exposing the case to extreme temperatures or harsh chemicals that may degrade its materials.
- **Training:** Ensure that you are trained in the proper use of the tourniquet and the RIGID TQ Case. Familiarize yourself with the deployment mechanism to ensure effective use in emergencies.

Instructions for Installation and Usage

1. Installation:

- Select a mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland drop leg holster.
- Use the provided mounting hardware to secure the RIGID TQ Case in the desired position.
- If using a BladeTech TEKLOK belt attachment, ensure it fits securely on belts up to 2.25 inches wide.
- Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.

2. Usage:

- To deploy the tourniquet, grasp the exposed portion of the tourniquet firmly.
- With a swift pull, the tourniquet will release from the case.
- Follow the manufacturer's instructions for applying the tourniquet effectively.

Disposal Instructions

- Dispose of the RIGID TQ Case in accordance with local regulations for plastic waste.
- If the case or tourniquet is damaged or expired, replace it immediately to ensure continued safety and effectiveness.

Contact Information for Further Support

For inquiries regarding safety or support, please refer to the contact information provided with your product packaging. Ensure you have the product details available for reference.

Thank you for prioritizing safety and preparedness with the Eleven 10 RIGID TQ Case for the CAT Gen 7.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk