RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, COYOTE

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10's products are designed, manufactured and assembled in the USA.

Attributes

- Name: RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, COYOTE
- Manufacturer: ELEVEN 10
- Product no.: 100043722
- Mfr. No.: E10-3011B-CYT
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668267

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RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions

Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This product is designed to securely hold your tourniquet, ensuring rapid access in emergency situations. Please read this safety instruction guide carefully to understand how to use and maintain your product safely.

General Safety Guidelines

- Ensure the RIGID TQ Case is used only for its intended purpose of holding the SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of damage or wear. Replace if necessary.
- Keep the product out of reach of children to prevent accidental misuse.
- Avoid exposing the case to extreme temperatures or harsh chemicals that could compromise its integrity.
- If you notice any unusual signs during use, such as difficulty in releasing the tourniquet, stop using the case immediately.

Specific Safety Precautions for Use

- Use the RIGID TQ Case only with compatible tourniquets (SOF®TT or SOF®TTW).
- Ensure that the case is securely mounted to your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- When deploying the tourniquet, ensure that your hands are positioned correctly to avoid any injuries.
- Be cautious while pulling the tourniquet from the case; ensure that you have a firm grip on the exposed portion to prevent slippage.
- If the case becomes wet, allow it to dry completely before use to maintain proper grip and function.

Instructions for Installation and Usage

1. Mounting the Case:

- Choose the desired mounting position (duty belt, MOLLE/PALS vest, or accessory position on a drop leg holster).
- If using the BladeTech TEKLOK belt attachment, ensure that it fits securely on a belt up to 2.25 inches wide.
- For MOLLE/PALS attachment, use the Tactical Tailor Malice Clip for secure mounting.
- If mounting to a Safariland® Drop Leg Shroud, use the appropriate mounting hardware (sold separately).

2. Securing the Tourniquet:

- Place the SOF®TT or SOF®TTW tourniquet inside the RIGID TQ Case.
- Ensure that the tourniquet is positioned in a way that allows for easy access and deployment.

3. Deploying the Tourniquet:

- To deploy, grasp the exposed portion of the tourniquet firmly.
- Pull swiftly to release the tourniquet from the case.
- Practice deployment in a safe environment to become familiar with the process.

4. Maintenance:

- Clean the case with a damp cloth as needed. Avoid using harsh chemicals.
- Check the mounting hardware regularly and use the included Vibratite to secure it in place once the desired configuration is obtained.

Disposal Instructions

- Dispose of the RIGID TQ Case responsibly. Check local regulations for proper disposal methods.
- If the case is damaged beyond repair, ensure it is disposed of in a manner that prevents accidental use.

Contact Information for Further Support

For any additional questions or concerns regarding the RIGID TQ Case, please refer to the manufacturer's website or customer support resources.

By following these safety instructions, you can ensure that your RIGID TQ Case is used effectively and safely. Thank you for your attention to these guidelines.

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