RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, RANGER GREEN

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10's products are designed, manufactured and assembled in the USA.

Attributes

- Name: RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, RANGER GREEN
- Manufacturer: ELEVEN 10
- Product no.: 100043723
- Mfr. No.: E10-3011B-RGR
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668298

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RIGID TQ Case for SOF®TT/SOF®TTW Safety Instruction Guide

Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This guide provides essential safety instructions and information to ensure the proper use and maintenance of your product. The RIGID TQ Case is designed for secure storage and rapid deployment of tourniquets, making it a crucial tool for first responders and individuals in emergency situations.

General Safety Guidelines

- Ensure the RIGID TQ Case is used only for its intended purpose: to hold and deploy the SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of damage or wear. If any damage is detected, discontinue use and replace the case.
- Keep the case away from extreme temperatures, moisture, and direct sunlight to maintain its structural integrity.
- Store the tourniquet in the case properly to ensure it is ready for immediate use.
- Familiarize yourself with the product and practice deploying the tourniquet in a safe environment before an emergency arises.

Specific Safety Precautions for Use

- Use the RIGID TQ Case only with compatible tourniquets (SOF®TT or SOF®TTW). Using incompatible products may result in failure to deploy the tourniquet effectively.
- Ensure the case is securely attached to your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster to prevent accidental detachment during use.
- Avoid using excessive force when pulling the tourniquet from the case to prevent injury or damage to the case.
- The case is designed for onehanded deployment; practice this technique to ensure proficiency, especially in highstress situations.
- Do not modify the case in any way. Alterations can compromise safety and functionality.

Instructions for Installation and Usage

1. Installation:

- Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- For Police Duty Belt: Use the BladeTech TEKLOK Belt Attachment, which fits up to a 2.25" belt.
- For MOLLE/PALS Attachment: Use the Tactical Tailor Malice Clip.
- Ensure the mounting hardware is securely fastened. Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.
- Refer to the installation instructions included with the mounting hardware for detailed steps.

2. Usage:

- To deploy the tourniquet, grasp the exposed portion of the tourniquet firmly.
- With a swift pull, remove the tourniquet from the case. The design allows for easy access even when hands are wet or bloody.
- Follow the manufacturer's instructions for applying the tourniquet in an emergency situation.

Disposal Instructions

- Dispose of the RIGID TQ Case responsibly. If the case is damaged or no longer needed, please recycle it according to local recycling regulations.
- Do not dispose of the case in regular household waste if it is made of materials that can be recycled.

Contact Information for Further Support

For any inquiries regarding safety, usage, or product support, please refer to the manufacturer's contact information provided on the product packaging or website.

Thank you for your attention to these safety guidelines. Your safety and the safety of others depend on your responsible use of this product. Stay safe!

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