PTAKS RUGGED REFILL POD - PTAKS RUGGED-REFILL POD W/COMBAT GAUZE COYOTE

Contents (1) CLEER Medical Trauma Dressing 4" (1) Phokus Research Frog Gauze OR QuikClot Combat Gauze LE (1) Nitrile Gloves (1) Phokus Research Group / Eleven 10 Frog Tape (1) North American Rescue Hyfin Vent Compact Chest Seal, Twin Pack Pouch Sold Seperate

Attributes

- Name: PTAKS RUGGED-REFILL POD W/COMBAT GAUZE COYOTE
- Manufacturer: ELEVEN 10
- Product no.: 100043754
- Mfr. No.: PTAKS-QC-CYT
- Size: 1 Person
- Delivery weight: 0.245kg
- Shipping height: 197mm
- Shipping width: 38mm
- Shipping length: 178mm
- UPC: 700615672356

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Safety Instruction Guide for PTAKS RUGGED REFILL POD

Introduction

Thank you for choosing the PTAKS RUGGED REFILL POD with Combat Gauze. This product is designed to provide essential first aid supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that all users are familiar with the contents and purpose of the PTAKS RUGGED REFILL POD.
- Store the product in a cool, dry place, away from direct sunlight and extreme temperatures.
- Regularly check the contents for expiration dates and signs of damage or tampering.
- Always use the product as intended and follow all instructions provided.
- Keep the product out of reach of children and vulnerable individuals when not in use.
- In case of any allergic reactions or adverse effects, discontinue use immediately and seek medical attention.

Specific Safety Precautions for Use

- Use the CLEER Medical Trauma Dressing only on clean, dry wounds.
- Ensure that Nitrile Gloves are worn when handling any bodily fluids to minimize contamination.
- The Phokus Research Frog Gauze or QuikClot Combat Gauze LE should be applied as directed to control bleeding.
- The North American Rescue Hyfin Vent Compact Chest Seal should be used according to the instructions to prevent air from entering the chest cavity.
- Dispose of any used materials safely and in accordance with local regulations.

Instructions for Installation and Usage

1. Preparation:

- Open the PTAKS RUGGED REFILL POD carefully.
- Review the contents to familiarize yourself with the items available for use.

2. Using the Trauma Dressing:

- Remove the dressing from its packaging.
- Place the dressing over the wound, ensuring complete coverage.
- Secure the dressing in place with tape or by applying pressure.

3. Applying Gauze:

- If bleeding is present, apply the Frog Gauze or QuikClot Combat Gauze LE directly to the wound.
- Apply firm pressure to control bleeding and hold the gauze in place.

4. Using the Chest Seal:

- Open the Hyfin Vent Compact Chest Seal package.
- Remove the adhesive backing and apply the seal over the wound on the chest.
- Ensure that the seal is airtight to prevent air from entering the chest cavity.

5. PostUse:

- After use, carefully dispose of all materials, including gloves, in a sealed bag.
- Wash your hands thoroughly after removing gloves.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations.
- Do not dispose of any items in regular household waste unless specified.
- Check for local hazardous waste disposal guidelines for items like gloves and dressings.

Contact Information for Further Support

For any inquiries regarding the safety or use of the PTAKS RUGGED REFILL POD, please contact your local distributor or the manufacturer directly. Ensure to have your product details ready for reference.

By following these safety instructions, you can ensure safe and effective use of the PTAKS RUGGED REFILL POD. Thank you for prioritizing safety and preparedness.

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