MBOK (MINI BLOW OUT KIT) KIT CONTENTS - MBOK KIT CONTENTS W/COMPRESSED GAUZE

Contents CLEER Medical Trauma Bandage (4" x 24" Flat Packed compression dressing) QuikClot Combat Gauze LE or Phokus Research Frog Gauze FOXSeal Occlusive Dressing (Contains 2 Chest Seals) Nitrile Gloves (1 pair) Phokus Research / Eleven 10® Frog Tape Pouch Sold Seperate

Attributes

- Name: MBOK KIT CONTENTS W/COMPRESSED GAUZE
- Manufacturer: ELEVEN 10
- Product no.: 100043758
- Mfr. No.: MBOK-GA
- Size: 1 Person
- Delivery weight: 0.227kg
- Shipping height: 51mm
- Shipping width: 178mm
- Shipping length: 203mm
- UPC: 700615674060

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MBOK (MINI BLOW OUT KIT) Safety Instruction Guide

Introduction

Thank you for choosing the MBOK (MINI BLOW OUT KIT). This kit is designed to provide essential first aid supplies for trauma situations. It is crucial to understand how to use this kit safely and effectively. This guide outlines general safety guidelines, specific precautions, installation and usage instructions, and disposal instructions to ensure safe and effective use of your MBOK kit.

General Safety Guidelines

- Ensure that all kit contents are intact and undamaged before use.
- Store the kit in a cool, dry place, away from direct sunlight and moisture.
- Regularly check the expiration dates of all medical supplies within the kit.
- Keep the kit out of reach of children to prevent misuse.
- In case of an emergency, call for professional medical assistance immediately.
- Report any unsafe products or incidents to the relevant authorities.

Specific Safety Precautions for Use

- Always wear nitrile gloves when handling blood or bodily fluids to reduce the risk of infection.
- Follow the instructions provided with each item in the kit for proper use.
- Use the CLEER Medical Trauma Bandage only on clean, dry wounds.
- When using QuikClot Combat Gauze or Phokus Research Frog Gauze, apply direct pressure to the wound to promote clotting.
- Ensure that FOXSeal Occlusive Dressing is applied correctly to prevent air from entering the chest cavity.
- Dispose of used medical supplies properly to prevent contamination or infection.

Instructions for Installation and Usage

Preparing the Kit for Use

1. **Inspect the Kit:** Before using the MBOK kit, inspect all items to ensure they are not expired or damaged.
2. **Gather Supplies:** Identify the supplies you will need based on the type of injury or emergency situation.

Using the Kit

- **Trauma Bandage:**

1. Remove the bandage from the packaging.
2. Place the bandage directly over the wound.
3. Secure it in place using the attached adhesive strips.

- **QuikClot Combat Gauze:**

1. Open the gauze package.
2. Pack the gauze into the wound.
3. Apply direct pressure for at least 3 minutes to promote clotting.

- **Phokus Research Frog Gauze:**

1. Open the gauze package.
2. Apply it directly to the wound, ensuring it is packed tightly.
3. Maintain pressure until bleeding stops.

- **FOXSeal Occlusive Dressing:**

1. Open the dressing package.
2. Remove the protective backing and place the seal over the chest wound.
3. Ensure it is sealed properly to prevent air from entering.

- **Nitrile Gloves:**

1. Put on gloves before handling any medical supplies.
2. Dispose of gloves after use to prevent contamination.

PostUse Instructions

1. **Clean Up:** Dispose of all used supplies in accordance with local regulations.
2. **Replenish the Kit:** Replace any used items to ensure the kit is ready for future use.

Disposal Instructions

- Dispose of used medical supplies in accordance with local health regulations.
- Do not throw away medical waste in regular trash.
- Check with local authorities for proper disposal methods for hazardous materials.

Contact Information for Further Support

For any inquiries or further assistance, please refer to the contact information provided with your product packaging or visit the manufacturer's website.

By following these guidelines, you can ensure the safe and effective use of your MBOK (MINI BLOW OUT KIT). Always prioritize safety and seek professional medical help when needed. Thank you for choosing our product, and we wish you the best in your preparedness efforts.

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