ETAK POUCH MEDICAL KIT CONTENTS - ETAK KIT CONTENTS W/COMBAT GAUZE

All the medical items to fill the ETAK Pouch or to replace expired products The Eleven 10 ETAK (Enhanced Trauma Aid Kit) is designed to hold critical medical items with the first responder and tactical officer in mind. Compact and durable, made with premium, duty-proven materials. For use with either the Fixed or Tear Away ETAK Pouches (sold separately). Kit Includes: (1) 4" Compact Medical Trauma Dressing (1) Compression Gauge or Combat Gauze (1) Nitrile Exam Gloves – Pair, Size L (1) Eleven 10® Frog Tape (1) NAR Hyfin Vent Compact Chest Seal, Twin Pack (1) 4" ChemLight Light Stick (1) 5.5" Trauma Shears

Attributes

- Name: ETAK KIT CONTENTS W/COMBAT GAUZE
- Manufacturer: ELEVEN 10
- Product no.: 100043775
- Mfr. No.: E10ETAKQC
- Size: 1 Person
- Delivery weight: 0.272kg
- Shipping height: 51mm
- Shipping width: 178mm
- Shipping length: 203mm
- UPC: 680306732951

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ETAK Pouch Medical Kit Safety Instructions

Introduction

Thank you for choosing the ETAK Pouch Medical Kit. This kit is designed to provide essential medical supplies for first responders and tactical officers. To ensure safe and effective use of this product, please read and follow the safety instructions outlined in this guide.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace expired items promptly.
- Keep the kit out of reach of children to prevent accidental use.
- Only use the kit for its intended purpose, as outlined in the instructions.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear protective gloves when handling medical supplies to reduce the risk of contamination.
- Follow the instructions for each item in the kit carefully to avoid misuse.
- If you are not trained in first aid or emergency response, seek assistance from a qualified individual when using the kit.
- In case of an emergency, prioritize calling emergency services before administering aid.

Instructions for Installation and Usage

1. Unpacking the Kit

- Carefully open the ETAK Pouch and inspect all contents.
- Ensure all items listed in the kit contents are present.

2. Using the Medical Items

- **4" Compact Medical Trauma Dressing**
 - Apply directly to the wound and secure in place.
- **Compression Gauge or Combat Gauze**
 - Use for severe bleeding; pack the gauze into the wound and apply pressure.
- **Nitrile Exam Gloves**
 - Put on a pair before handling any medical supplies.
- **Eleven 10@ Frog Tape**
 - Use to secure dressings or bandages in place.
- **NAR Hyfin Vent Compact Chest Seal**
 - Apply to any open chest wound to prevent air from entering the chest cavity.
- **4" ChemLight Light Stick**
 - Activate in lowlight situations for visibility and safety.
- **5.5" Trauma Shears**
 - Use to cut clothing or bandages as needed.

3. PostUse Care

- Dispose of used medical items according to local regulations.
- Clean the kit and restock any used or expired items as soon as possible.

Disposal Instructions

- Dispose of any expired or used medical items in accordance with local waste management regulations.
- Do not dispose of medical items in regular household waste unless specified by local guidelines.
- Recycle the pouch if it is made from recyclable materials; check local recycling rules.

Contact Information for Further Support

For any safety inquiries or further assistance regarding the ETAK Pouch Medical Kit, please refer to the contact information provided on the packaging or the manufacturer's website.

By following these safety instructions, you can ensure the effective use of your ETAK Pouch Medical Kit while adhering to safety regulations. Thank you for your attention to these important guidelines.

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