MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MOUNT BLACK

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MOUNT BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050370
- Mfr. No.: BT-TKN-MTKN-ESS-BK
- Size: 1 Person
- Delivery weight: 0.318kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650107

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Micro Trauma Kit NOW! Essentials Supplies Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This compact and efficient trauma kit is designed for everyday carry, providing essential lifesaving medical supplies. This guide outlines important safety information to ensure safe usage and compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Always inspect the kit before use to ensure all components are intact and undamaged.
- Keep the kit out of reach of children unless supervised by an adult.
- Be aware of the specific uses of each item in the kit and follow instructions for use.
- In case of an emergency, use the kit as intended and seek professional medical assistance immediately.
- Regularly check for updates on product recalls and safety notices through the EU Safety Gate platform.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Ensure gloves are intact before use to avoid contamination.
 - Dispose of used gloves properly to prevent infection.
- Quickclot Bleeding Control Dressing:
 - Apply directly to the wound; do not use on nonbleeding areas.
 - Seek medical help for severe injuries.
- Mini Compression Dressing:
 - Use as directed to control bleeding; do not apply too tightly.
- Hypothermia Blanket:
 - Use to maintain body temperature; avoid direct contact with fire or heat sources.
- Flat Fold Tape:
 - Ensure the area is clean before applying tape; avoid sensitive skin.

Instructions for Installation and Usage

- **Belt Mounting:**
 1. Locate the MOLLE or belt mounting system on the pouch.
 2. Attach the pouch securely to your belt or MOLLE system, ensuring it is stable and accessible.
- **Accessing Supplies:**
 1. Use the BLIP pull tabs to access supplies quickly with one hand or a single finger.
 2. Practice accessing the kit to ensure efficiency during emergencies.
- **Maintaining the Kit:**
 - Regularly check the contents of the kit to ensure supplies are stocked and within expiration dates.
 - Replace any used or expired items promptly.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not throw away the kit itself unless it is damaged beyond repair.
- Recycle the pouch if possible, following local recycling guidelines.

Contact Information for Further Support

For any safety inquiries or concerns regarding the Micro Trauma Kit NOW!, please refer to the contact information provided with your product packaging or visit the official website for updates and support.

Thank you for prioritizing safety and preparedness with the Micro Trauma Kit NOW! Always remember to be aware of your surroundings and seek professional help in emergencies.

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