MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-MULTICAM

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100988
- Mfr. No.: TKN-NANO-ESS-MC
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073652835

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact trauma pouch is designed to provide essential medical supplies for emergency situations. It is critical to understand how to use this kit safely and effectively, as it is intended to help manage traumatic injuries, particularly bleeding. Please read this safety instruction guide carefully to ensure proper usage and compliance with safety regulations.

General Safety Guidelines

- Always check the contents of the kit before use to ensure all items are present and in good condition.
- Store the kit in a cool, dry place away from direct sunlight to maintain the integrity of the medical supplies.
- Keep the kit out of reach of children to prevent accidental use or ingestion of any items.
- Be aware of potential allergies to materials used in the kit, such as latex in gloves.
- In case of any adverse reactions or injuries, seek medical assistance immediately.
- Familiarize yourself with basic first aid principles to maximize the effectiveness of the kit.

Specific Safety Precautions for Use

- **Trauma Gloves:**
 - Always wear gloves when treating wounds to prevent infection and protect yourself from bloodborne pathogens.
 - Dispose of gloves safely after use to avoid contamination.
- **QuikClot Products:**
 - Follow the instructions provided with the QuikClot products carefully to ensure effective bleeding control.
 - Do not use QuikClot products on nonbleeding wounds or for purposes other than indicated.
- **Pressure Device:**
 - Ensure the pressure device is applied correctly to avoid further injury.
 - Monitor the injured area for signs of excessive swelling or complications.

Instructions for Installation and Usage

1. Preparation:

- Open the Micro Trauma Kit NOW! Nano and review the contents.
- Identify the items you will need based on the nature of the injury.

2. Using Trauma Gloves:

- Roll out one pair of trauma gloves and put them on before handling any injured person or material.

3. Applying QuikClot Dressing:

- If bleeding is present, unpack the QuikClot Bleeding Control Dressing.
- Apply the dressing directly to the wound, ensuring it covers the entire area.
- Hold the dressing in place with gentle pressure for several minutes.

4. Using the Pressure Device:

- If necessary, use the pressure device to maintain consistent pressure on the dressing.
- Ensure the device is secured but not so tight that it cuts off circulation.

5. Aftercare:

- Once bleeding is controlled, seek further medical assistance.
- Remove gloves and dispose of them properly.

Disposal Instructions

- Dispose of used medical supplies, including gloves and dressings, in a sealed plastic bag to prevent contamination.
- Follow local regulations for the disposal of medical waste.
- Do not reuse any items from the kit that have been used or are damaged.

Contact Information for Further Support

For any safety inquiries or concerns regarding the Micro Trauma Kit NOW! Nano, please reach out to the designated EU contact point provided with your product packaging.

Thank you for taking the time to read this safety instruction guide. Your safety and the safety of others depend on the proper use of this kit. Always be prepared and informed when it comes to emergency situations.

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