MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-RANGER GREEN

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100989
- Mfr. No.: TKN-NANO-ESS-RG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 25mm
- Shipping width: 178mm
- Shipping length: 279mm
- UPC: 810073652842

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MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! This compact trauma kit is designed to provide essential medical supplies for emergency situations. It is crucial to use this kit safely and effectively to ensure the best outcomes in case of trauma. This guide outlines important safety instructions in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Always read the instructions carefully before using the kit.
- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure that all items are intact and within their expiration dates.
- Keep the kit out of reach of children and vulnerable individuals to prevent accidental use.
- If you experience any allergic reactions to the materials in the kit, seek medical attention immediately.
- Report any unsafe products or accidents related to this kit to the appropriate authorities.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Use gloves when handling blood or bodily fluids to minimize the risk of infection.
 - Dispose of used gloves in a safe manner to prevent contamination.
- QuikClot Bleeding Control Dressing:
 - Apply the dressing directly to the wound and apply pressure as needed.
 - Do not remove the dressing once applied, as this may disrupt clotting.
- Pressure Device:
 - Ensure the pressure device is applied correctly to control bleeding effectively.
 - Monitor the patient for any signs of shock or distress.

Instructions for Installation and Usage

- **How to Use the Kit:**
 1. **Assess the Situation:**
 - Ensure your safety and the safety of the injured person.
 - Call for emergency medical assistance if necessary.
 2. **Prepare the Area:**
 - Find a clean and safe area to treat the injury.
 - Gather all necessary supplies from the kit.
 3. **Wear Trauma Gloves:**
 - Put on the trauma gloves before handling any medical supplies or the injured person.
 4. **Control Bleeding:**
 - Use the QuikClot Bleeding Control Dressing:
 - Open the dressing and apply it directly to the wound.
 - Apply firm pressure for at least 3 minutes.
 5. **Apply Pressure Device:**
 - If bleeding persists, use the pressure device to maintain pressure on the dressing.
 - Ensure it is secured properly but not too tight to cut off circulation.
 6. **Monitor Patient:**
 - Keep the patient calm and monitor their condition until help arrives.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local regulations.
- Ensure that all contaminated materials, such as gloves and dressings, are disposed of in biohazard waste containers if available.
- Regularly check the kit for expired items and dispose of them safely.

Contact Information for Further Support

For any inquiries regarding the product, please refer to the contact information provided by the retailer or manufacturer. They will be able to assist you with any safety concerns or product support.

By following these safety instructions, you can ensure the effective use of the MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! while adhering to the EU General Product Safety Regulation. Always prioritize safety and seek professional medical assistance when necessary.

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