MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-WOLF GRAY

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100990
- Mfr. No.: TKN-NANO-ESS-WF
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 38mm
- Shipping width: 89mm
- Shipping length: 171mm
- UPC: 810073652859

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MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! Nano. This compact trauma pouch is designed to provide essential medical supplies for emergency situations, specifically to control bleeding, which is a leading cause of traumatic death. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the kit is kept in a dry and cool place to maintain the integrity of the medical supplies.
- Regularly check the contents of the kit to ensure that all items are present and not expired.
- Keep the kit out of reach of children unless under adult supervision.
- In case of an emergency, always call for professional medical assistance after administering first aid.
- Report any unsafe products or incidents related to this kit to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to minimize the risk of infection.
- Use the QuikClot bleeding control dressing as directed to effectively control bleeding.
- Do not use any items in the kit if they are damaged or expired.
- If you are unsure about how to use any of the supplies, seek guidance from a qualified medical professional.
- Be cautious when applying pressure with the pressure device; excessive pressure can cause further injury.

Instructions for Installation and Usage

- **Preparation:**
 - Open the kit and lay out all components on a clean surface.
 - Familiarize yourself with each item and its intended use.
- **Using the Trauma Gloves:**
 - Roll out the gloves from their packaging.
 - Ensure hands are clean before putting on the gloves.
 - Dispose of gloves properly after use.
- **Applying the QuikClot Bleeding Control Dressing:**
 - Identify the source of bleeding.
 - Remove the dressing from its packaging.
 - Place the dressing directly on the wound and apply firm pressure for at least 3 minutes.
 - If bleeding continues, apply a second dressing on top without removing the first.
- **Using the Pressure Device:**
 - Position the pressure device over the dressing.
 - Secure it in place according to the product instructions.
 - Monitor the wound and adjust pressure as necessary, ensuring it is not too tight.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not dispose of medical waste in regular household trash.
- Follow guidelines for hazardous waste disposal if applicable.

Contact Information for Further Support

For safety inquiries or additional support related to the MICRO TRAUMA KIT NOW! Nano, please refer to the product packaging for contact details of the manufacturer or authorized distributor.

Thank you for prioritizing safety and being prepared for emergencies. Your MICRO TRAUMA KIT NOW! Nano is an important tool in providing immediate care when it is needed most.

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