MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-BLACK

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100991
- Mfr. No.: TKN-NANO-PRO-BK
- Size: 1 Person
- Delivery weight: 0.227kg
- Shipping height: 121mm
- Shipping width: 83mm
- Shipping length: 44mm
- UPC: 810073652866

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact first aid kit is designed to provide essential medical supplies for managing traumatic injuries, particularly severe bleeding. Please read this safety instruction guide carefully to ensure the safe and effective use of this product.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items immediately.
- Keep the kit out of reach of children to prevent accidental ingestion or misuse.
- Use the kit only for its intended purpose as outlined in this guide.
- In case of severe injuries or emergencies, seek professional medical assistance immediately.

Specific Safety Precautions for Use

- Always wear Trauma Gloves when handling any medical supplies to minimize the risk of contamination.
- Use Quikclot products according to the instructions provided on the packaging to ensure effective bleeding control.
- Be aware of potential allergic reactions to materials in the kit, such as latex in gloves. If you have allergies, consider using nonlatex alternatives.
- If you are unsure how to use any item in the kit, seek guidance from a qualified medical professional before attempting to use it.
- Do not use the kit for treating nontraumatic injuries or conditions that require specialized medical attention.

Instructions for Installation and Usage

1. Opening the Kit:

- Open the pouch carefully to avoid damaging any contents.
- Lay out all items on a clean surface for easy access.

2. Using Trauma Gloves:

- Roll out one pair of Trauma Gloves from the pouch.
- Put on the gloves before handling any medical supplies.

3. Applying Quikclot Bleeding Control Dressing:

- Identify the source of bleeding.
- Open the Quikclot dressing package.
- Apply the dressing directly to the wound, ensuring it covers the entire area.
- Apply pressure for at least 3 minutes to promote clotting.

4. Using Pressure Device:

- If bleeding does not stop with the dressing, use the Pressure Device as instructed on its packaging.
- Ensure it is applied firmly to maintain pressure on the wound.

5. PostUse Care:

- Dispose of used gloves and medical supplies in accordance with local regulations.
- Wash your hands thoroughly after removing the gloves.

Disposal Instructions

- Dispose of all used items in accordance with local waste disposal regulations.
- Do not dispose of the kit in regular household trash if it contains hazardous materials.
- Check local guidelines for disposing of medical waste to ensure safe and compliant disposal.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW! Nano, including safety concerns or product recalls, please contact the appropriate regulatory authority in your country. Ensure that you have the product details at hand for a swift response.

By following these guidelines, you can ensure the safe and effective use of your Micro Trauma Kit NOW! Nano. Thank you for prioritizing safety and preparedness.

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