MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-COYOTE BROWN

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-COYOTE BROWN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100992
- Mfr. No.: TKN-NANO-PRO-CB
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 38mm
- Shipping width: 83mm
- Shipping length: 178mm
- UPC: 810073652873

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MICRO TRAUMA KIT NOW! NANO USER MANUAL

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! NANO. This compact trauma kit is designed to provide essential medical supplies for emergency situations. It is vital for stopping traumatic bleeding, which is the leading cause of death in such scenarios. Please read this manual carefully to understand how to use the kit safely and effectively.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and damage.
- Keep the kit out of reach of children unless supervised by an adult.
- Use the kit only for its intended purpose and according to the instructions provided.
- In case of an emergency, seek professional medical assistance immediately after using the kit.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling any medical supplies to maintain hygiene.
- Ensure that the QuikClot products are used in accordance with the instructions provided in the packaging.
- Do not use any item in the kit if the packaging is damaged or if the item appears to be compromised.
- If you are uncertain about how to use any item, seek professional guidance before attempting to use it.

Instructions for Installation and Usage

1. Opening the Kit

- Unzip or open the pouch carefully to avoid damaging the contents.
- Familiarize yourself with the layout of the kit and the location of each item.

2. Using Trauma Gloves

- Roll the gloves out of their packaging.
- Ensure your hands are clean before putting on the gloves.
- Dispose of gloves properly after use to maintain hygiene.

3. Applying QuikClot Products

- For the QuikClot Bleeding Control Dressing:
 - Remove the dressing from its packaging.
 - Place the dressing directly on the wound and apply firm pressure.
 - Secure the dressing in place with additional bandaging if necessary.
- For the QuikClot Combat Gauze:
 - Open the gauze and apply it directly to the bleeding site.
 - Pack the gauze into the wound if necessary and apply pressure.
 - Maintain pressure for several minutes until bleeding stops.

4. Using the Pressure Device

- Follow the instructions provided with the pressure device for proper application.
- Ensure that the device is securely fastened to maintain pressure on the wound.

5. PostUse Care

- After using any items, remove gloves and dispose of them properly.
- Clean the area around the wound if possible and seek further medical assistance.

Disposal Instructions

- Dispose of used items in accordance with local regulations regarding medical waste.
- Ensure that all used materials, including gloves and dressings, are sealed in a plastic bag before disposal to prevent contamination.
- Check local guidelines for disposal of expired or unused items from the kit.

Contact Information for Further Support

For any inquiries regarding safety or product usage, please reach out to the designated EUbased contact point as provided on the packaging or product website.

By following these guidelines and instructions, you can ensure the safe and effective use of the MICRO TRAUMA KIT NOW! NANO. Remember, while this kit is designed to assist in emergencies, it is not a substitute for professional medical care. Always seek help from qualified medical personnel when needed.

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