MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-MULTICAM

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100993
- Mfr. No.: TKN-NANO-PRO-MC
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073652880

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact medical kit is designed to provide essential first aid supplies for emergency situations. It is important to understand how to use this product safely and effectively. This guide outlines the safety instructions and precautions necessary for the proper use of this kit in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure that all contents of the kit are intact and not expired before use.
- Store the kit in a cool, dry place away from direct sunlight.
- Keep the kit out of reach of children unless under adult supervision.
- Regularly check the kit for any missing or damaged items and replace them as necessary.
- If you experience an allergic reaction to any of the materials in the kit, discontinue use immediately and seek medical attention.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to prevent infection.
- Use QuikClot products only as directed on the packaging. Do not apply to large or deep wounds without medical supervision.
- Ensure that the pressure device is applied correctly to avoid excessive pressure that may cause further injury.
- If you are unsure about how to use any item in the kit, seek guidance from a qualified medical professional.
- Do not use the kit for purposes other than first aid treatment.

Instructions for Installation and Usage

1. Preparation:

- Open the pouch and lay out all items to ensure you have everything needed for the situation.
- Familiarize yourself with each item and its intended use.

2. Applying QuikClot Bleeding Control Dressing:

- Use gloves before handling any dressings.
- Apply the dressing directly to the wound, ensuring it covers the entire area.
- Hold in place with gentle pressure until bleeding stops.

3. Using the Pressure Device:

- Position the device over the dressing.
- Secure it in place, ensuring it is snug but not overly tight.
- Monitor the wound for any signs of excessive bleeding or discomfort.

4. Disposing of Used Items:

- Place used gloves, dressings, and any other contaminated materials in a sealed plastic bag.
- Dispose of the bag in accordance with local regulations for biohazard waste.

5. PostUse Care:

- After use, clean your hands thoroughly.
- Restock any items used from the kit to ensure it is ready for future emergencies.

Disposal Instructions

- Dispose of the Micro Trauma Kit NOW! and its contents in accordance with local regulations.
- Do not dispose of the kit in regular household waste if it contains biohazard materials.
- Check for any specific disposal instructions on individual items within the kit.

Contact Information for Further Support

For any inquiries or concerns regarding the Micro Trauma Kit NOW!, please refer to your local distributor or check the packaging for additional contact details.

By following these safety guidelines and instructions, you can ensure the effective and safe use of the Micro Trauma Kit NOW! Nano. Always prioritize safety and seek professional medical assistance when necessary.

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