MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-WOLF GRAY

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-PRO SUPPLIES-WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100995
- Mfr. No.: TKN-NANO-PRO-WF
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 38mm
- Shipping width: 83mm
- Shipping length: 165mm
- UPC: 810073652903

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano, the smallest TKN medical kit designed to provide essential medical supplies for managing traumatic injuries. This guide outlines important safety information and instructions for the safe use and disposal of this product, ensuring compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure the product is used only for its intended purpose: to provide first aid for traumatic injuries.
- Keep the kit out of reach of children to prevent accidental use or ingestion of contents.
- Regularly inspect the kit for any damaged or expired components before use.
- Store the kit in a cool, dry place to maintain the integrity of the medical supplies.
- In case of an emergency, call for professional medical assistance immediately.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Use gloves to protect yourself and the injured person from infection.
 - Ensure gloves are properly fitted and free of tears before use.
- QuikClot Bleeding Control Dressing:
 - Apply the dressing directly to the wound to control bleeding.
 - Do not remove the dressing once applied; it is designed to be left in place.
- Pressure Device:
 - Use the pressure device to apply consistent pressure to the wound.
 - Ensure the device is securely fastened to prevent it from slipping during use.

Instructions for Installation and Usage

1. Preparation:

- Before using the kit, ensure you are in a safe environment.
- Wash your hands thoroughly or use hand sanitizer if available.

2. Using Trauma Gloves:

- Open the glove package carefully.
- Put on the gloves before attending to the injured person.

3. Applying QuikClot Dressing:

- Open the QuikClot package.
- Place the dressing directly onto the bleeding wound.
- Hold the dressing in place for at least 3 minutes to allow clotting.

4. Using the Pressure Device:

- Secure the pressure device over the QuikClot dressing.
- Adjust the device to apply firm pressure without cutting off circulation.
- Monitor the injured person for any changes in condition.

5. After Use:

- Dispose of used gloves and any other waste materials in accordance with local regulations.
- Clean any reusable items according to the manufacturer's instructions.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local regulations for hazardous waste.
- Do not reuse single-use items.
- Store the kit in a safe location until the next use, ensuring all contents are intact.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW! Nano, please refer to the manufacturer's website or contact customer support for assistance. Remember to report any unsafe products or incidents to the appropriate authorities.

Conclusion

By following these safety instructions, you can ensure the effective and safe use of your Micro Trauma Kit NOW! Nano. Always prioritize safety and seek professional medical assistance when needed. Thank you for your commitment to safety and preparedness.

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