MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-BLACK

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100996
- Mfr. No.: TKN-NANO-ADV-BK
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073652910

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact medical kit is designed to provide essential supplies for managing traumatic injuries, particularly those involving bleeding. This guide will help you understand how to safely use and maintain your kit, ensuring you are prepared in case of an emergency.

General Safety Guidelines

- Always read and understand the instructions before using the kit.
- Keep the kit out of reach of children unless under adult supervision.
- Regularly check the contents of the kit to ensure items are not expired or damaged.
- Store the kit in a cool, dry place to maintain the integrity of the supplies.
- In case of an emergency, seek professional medical help immediately after using the kit.

Specific Safety Precautions for Use

- **Trauma Gloves:**
 - Use gloves to protect yourself from blood and bodily fluids.
 - Ensure gloves are intact and free from tears before use.
- **QuikClot Dressing:**
 - Apply the dressing directly to the bleeding wound.
 - Do not remove the dressing once applied; it is designed to be left in place until medical help arrives.
- **Pressure Device:**
 - Use the pressure device to apply direct pressure to the wound.
 - Ensure it is securely fastened but not so tight that it cuts off circulation.

Instructions for Installation and Usage

1. Preparation:

- Ensure you are in a safe environment before administering first aid.
- Wash your hands thoroughly, or use hand sanitizer if soap and water are unavailable.

2. Using the Kit:

- Open the kit and identify the items you need based on the injury.
- Put on the trauma gloves before handling any medical supplies.
- For bleeding control:
 - Take the QuikClot dressing and apply it directly to the wound.
 - If the bleeding continues, use the pressure device for additional pressure.
- After applying the necessary supplies, monitor the injured person for any changes in condition until professional help arrives.

3. PostUse:

- Dispose of any used gloves and medical supplies in accordance with local regulations.
- Replenish any used items in the kit as soon as possible to ensure it is ready for future use.

Disposal Instructions

- Dispose of used gloves and medical supplies in a sealed plastic bag to prevent contamination.
- Follow local regulations for the disposal of medical waste.
- Regularly check the kit for expired items and dispose of them properly.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW! Nano, please refer to the manufacturer's contact information included with the product packaging. It is essential to have access to support for safety inquiries and product recalls.

By following these safety instructions, you can ensure that you are using the Micro Trauma Kit NOW! Nano effectively and safely. Always remember that this kit is a first aid tool and should not replace professional medical assistance.

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