MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-COYOTE BROWN

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-COYOTE BROWN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100997
- Mfr. No.: TKN-NANO-ADV-CB
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 38mm
- Shipping width: 95mm
- Shipping length: 171mm
- UPC: 810073652927

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact first aid kit is designed to provide essential medical supplies for trauma situations, particularly to control bleeding, which is the leading cause of traumatic death. Please read this safety instruction guide carefully to ensure the safe and effective use of this product.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Keep the kit out of reach of children to prevent accidental use.
- Regularly check the contents of the kit to ensure all items are intact and within their expiration dates.
- Do not use any items that appear damaged or compromised.
- Familiarize yourself with the contents and their intended use before an emergency arises.

Specific Safety Precautions for Use

- Always wear trauma gloves when providing first aid to prevent infection and contamination.
- Use QuikClot products strictly according to the instructions provided with each item.
- If you are unsure about how to use any item in the kit, seek professional medical assistance immediately.
- In case of severe bleeding, apply pressure directly to the wound and seek emergency medical help as soon as possible.
- Keep the kit accessible, but ensure it is stored securely to prevent unauthorized access.

Instructions for Installation and Usage

1. Unpacking the Kit:

- Carefully open the kit and inspect all items for damage.
- Ensure that all essential components are present: Trauma Gloves, Quikclot Bleeding Control Dressing, and Pressure Device.

2. Using the Trauma Gloves:

- Remove the gloves from their packaging.
- Put on the gloves before attending to any injuries to maintain hygiene.
- Dispose of the gloves properly after use to prevent contamination.

3. Applying QuikClot Bleeding Control Dressing:

- Open the dressing package carefully.
- Place the dressing directly over the bleeding wound.
- Apply firm pressure for at least 3 minutes while monitoring the bleeding.
- If bleeding continues, do not remove the dressing; add additional dressings on top.

4. Using the Pressure Device:

- Follow the specific instructions provided with the pressure device.
- Apply the device to the wound as directed, ensuring it is secure but not overly tight.

5. Aftercare:

- Once the bleeding is controlled, seek professional medical assistance.
- Monitor the injured person for signs of shock or other complications.

Disposal Instructions

- Dispose of used items, such as gloves and dressings, in accordance with local regulations for hazardous waste.
- Ensure that the kit is replenished after use to maintain readiness for future emergencies.
- Check expiration dates regularly and replace expired items promptly.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW! Nano, please reach out to your local distributor or visit the official website for more information.

Important Note: This guide adheres to the EU General Product Safety Regulation (GPSR) and aims to provide consumers with essential safety information regarding the use of the Micro Trauma Kit NOW! Nano. Always prioritize safety and seek professional help in emergencies.

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