MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-RANGER GREEN

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ADVANCED SUPPLIES-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100999
- Mfr. No.: TKN-NANO-ADV-RG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 25mm
- Shipping width: 178mm
- Shipping length: 279mm
- UPC: 810073652941

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MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! Safety Instructions

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! Nano. This compact trauma kit is designed to provide essential medical supplies for emergency situations. It is crucial to understand how to use this product safely and effectively. This guide outlines important safety information, usage instructions, and disposal guidelines to ensure your safety and the safety of others.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the kit for expired items and replace them as necessary.
- Keep the kit out of reach of children unless under adult supervision.
- Always read the instructions before use to familiarize yourself with the contents and their applications.
- In case of a medical emergency, seek professional medical assistance immediately after using the kit.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Always wear gloves when handling blood or bodily fluids to prevent contamination.
 - Dispose of used gloves properly to avoid crosscontamination.
- QuikClot Bleeding Control Dressing:
 - Apply directly to the wound with pressure to control bleeding.
 - Do not remove the dressing once applied; seek professional medical help.
- Pressure Device:
 - Use as directed to apply pressure to bleeding wounds.
 - Ensure the device is securely fastened but not overly tight to avoid restricting blood flow.

Instructions for Installation and Usage

1. Preparation:

- Open the Micro Trauma Kit NOW! Nano and lay out all items.
- Ensure you have a clear understanding of each item's purpose.

2. Using Trauma Gloves:

- Open the glove package carefully.
- Put on the gloves before handling any medical supplies or treating wounds.

3. Applying QuikClot Dressing:

- Identify the source of bleeding.
- Place the QuikClot dressing directly onto the wound.
- Apply firm pressure for at least 3 minutes to allow clotting.

4. Using the Pressure Device:

- Position the pressure device over the dressing.
- Secure it in place, ensuring it applies pressure without cutting off circulation.
- Monitor the area for signs of continued bleeding or swelling.

5. After Use:

- Carefully remove and dispose of gloves and other used materials.
- Clean your hands thoroughly after assisting with any medical treatment.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local regulations for hazardous waste.
- Do not dispose of items in regular household waste unless specified by local guidelines.
- Keep the kit container for future use; it can be restocked with new supplies after use.

Contact Information for Further Support

For any inquiries regarding the MICRO TRAUMA KIT NOW! Nano, including safety concerns or product recalls, please refer to the appropriate channels for assistance. Always ensure to stay updated with the latest safety information through official resources.

Conclusion

The MICRO TRAUMA KIT NOW! Nano is an essential tool for emergency preparedness. By following these safety instructions and guidelines, you can ensure safe and effective use of your trauma kit. Remember, in any emergency situation, professional medical help should always be sought after initial first aid measures. Stay safe and be prepared!

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