

DELTA BAG - SPIRITUS SYSTEMS DELTA BAG MULTICAM ARID

[Emergency & Survival Gear](#) > [Backpacks & Bags](#) > [Backpacks](#)

The Delta Bag is a compact, assault medic bag that is designed to be worn over a plate carrier and hold a wide variety of medical supplies and equipment. The front of the bag is divided into three uniformly sized compartments with diagonal zipper access that enables the end user to quickly access the entire pocket without having to reposition any of its contents. The unique configuration also aids in passive retention of gear when the compartment is opened. The main compartment of the Delta Bag is a clamshell design with two full fields of loop that can accept any hookbacked accessories (including the entire line of Spiritus Systems utility inserts, and other 3rd party products). The bag can be transported with the clamshell in the open position from location to location using the integrated briefcase handles, and provides an intuitive workstation when emergency medical care is being administered. The back side of the pack features two slim compartments for keeping flat-pack items secure and damage free, and keeps the pack close to the body when worn. The included straps for the delta bag are removeable and replaceable with any 3rd party straps that attach with 1" webbing. The bag features four sewn-in metal tri-glide mounting points for strap integration. The Delta bag features mesh and metal eyelets in key locations for full water drainage capability, and is designed to remain effective and accessible in austere environments.



Attributes

- Name: [SPIRITUS SYSTEMS DELTA BAG MULTICAM ARID](#)
- Manufacturer: [SPIRITUS SYSTEMS](#)
- Product no.: 430104964
- Mfr. No.: BU-DB-MA
- Colour: Multi-Cam Arid
- Weight: 1.261 lbs
- Delivery weight: 0.587kg
- Shipping height: 76mm
- Shipping width: 457mm
- Shipping length: 305mm

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Delta Bag Spiritus Systems Safety Instruction Guide](#)
- [About Us](#)

Delta Bag Spiritus Systems Safety Instruction Guide

Introduction

Thank you for choosing the Delta Bag by Spiritus Systems. This compact assault medic bag is designed for efficient use in various environments. To ensure your safety and the optimal performance of the product, please read and follow these safety instructions carefully.

General Safety Guidelines

- Always inspect the Delta Bag before use for any signs of damage or wear.
- Ensure that all medical supplies and equipment stored in the bag are within their expiration dates and properly secured.
- Use the bag only for its intended purpose as a medical supply carrier.
- Keep the bag out of reach of children when not in use.
- Do not overload the bag beyond its intended capacity to avoid compromising its integrity.
- Be aware of your surroundings when using the bag in emergency situations.

Specific Safety Precautions for Use

- **Zipper Access:** Always ensure that zippers are fully closed when the bag is not in use to prevent accidental opening.
- **Weight Distribution:** Distribute weight evenly within the compartments to maintain balance and comfort when worn.
- **Strap Usage:** When using the included straps, ensure they are securely attached to the mounting points to prevent accidental detachment.
- **Water Drainage:** Utilize the mesh and metal eyelets for water drainage, especially in wet conditions, to keep contents dry and maintain bag integrity.
- **Transporting the Bag:** When transporting the bag in an open clamshell position, ensure that it is securely held to prevent accidental spillage of contents.

Instructions for Installation and Usage

1. Strap Attachment:

- Locate the four sewn-in metal triglide mounting points on the bag.
- Attach the included straps or any compatible 1" webbing straps securely to the mounting points.

2. Packing the Bag:

- Open the bag using the diagonal zippers.
- Organize medical supplies within the three compartments for easy access.
- Use hookbacked accessories to secure additional items in the main compartment.

3. Wearing the Bag:

- Adjust the straps for a comfortable fit.
- Ensure the bag is close to your body to prevent shifting during movement.

4. Using the Bag in Emergencies:

- Open the bag fully to access all compartments.

- Use the bag as a workstation by placing it on a stable surface if needed.

Disposal Instructions

- Dispose of the Delta Bag in accordance with local regulations for textile and plastic waste.
- If the bag becomes damaged beyond repair, consider recycling components where possible.

Contact Information for Further Support

For any safety inquiries or additional support, please refer to the manufacturer's contact details provided in the product packaging or visit the official website.

Please keep this safety instruction guide for future reference. Your safety and the effective use of the Delta Bag are our top priorities.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk