

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, RANGER GREEN

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/TEK-LOK, RANGER GREEN](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043713
- Mfr. No.: E10-3021B-RGR
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670185

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Rigid TQ Case for CAT with TEKLOK Safety Instruction Guide](#)
- [About Us](#)

Rigid TQ Case for CAT with TEKLOK Safety Instruction Guide

Introduction

Thank you for choosing the Eleven 10® RIGID TQ Case® for the CAT® Gen 7. This hardsided tourniquet case is designed for easy access and rapid deployment of your tourniquet in emergency situations. Please read this safety instruction guide carefully to ensure safe use and proper handling of the product.

General Safety Guidelines

- Ensure that the RIGID TQ Case® is used only for the purpose it is designed for, which is to securely hold and provide quick access to the CAT® Tourniquet.
- Regularly inspect the case for any signs of wear or damage. Do not use the case if it is compromised.
- Keep the case and its contents away from children and vulnerable individuals, as improper use may lead to accidents.
- Follow all local laws and regulations regarding the use of medical equipment and personal protective gear.

Specific Safety Precautions for Use

- Always ensure that the tourniquet is properly secured within the RIGID TQ Case® before use.
- When deploying the tourniquet, do so with one hand using gross motor skills to ensure quick access, especially in highpressure situations.
- Avoid exposing the case to extreme temperatures or harsh chemicals, as this may affect its integrity and functionality.
- If the case is mounted on a duty belt or vest, ensure that it is securely attached to prevent accidental removal or loss during use.

Instructions for Installation and Usage

Installation

- Select a mounting position that is easily accessible for rapid deployment (vertically, horizontally, canted, or inverted).
- Use the provided hardware to secure the RIGID TQ Case® to your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained. Instructions for application are included with the product.

Usage

- To deploy the tourniquet:
 - Grasp the exposed portion of the tourniquet.
 - With a swift pull, release the tourniquet from the case.
- After use, ensure that the tourniquet is reset and returned to the case for future use.
- Regularly practice deploying the tourniquet to maintain proficiency, especially in emergency scenarios.

Disposal Instructions

- Dispose of the RIGID TQ Case® and any used tourniquets in accordance with local regulations for medical waste.
- Do not incinerate the case, as this may produce harmful fumes.

- If the case is no longer functional or has reached the end of its usable life, please recycle it according to local recycling guidelines where applicable.

Contact Information for Further Support

For any inquiries regarding safety, usage, or product support, please consult the contact information provided with your product packaging or visit the manufacturer's website.

Thank you for your attention to these safety instructions. Your safety and the effective use of the RIGID TQ Case® are our top priorities.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk