

RIGID TQ CASE FOR C-A-T WITH TEK-LOK - ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID CAT GEN 7 CASE W/MALICE CLIP, MULTICAM](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043720
- Mfr. No.: E10-3021M-MTC
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615672035

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Rigid TQ Case for CAT with TekLok Safety Instructions](#)
- [About Us](#)

Rigid TQ Case for CAT with TekLok Safety Instructions

Introduction

Thank you for choosing the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. This case is designed to provide quick and reliable access to your tourniquet in emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the Rigid TQ Case is used only for its intended purpose.
- Regularly inspect the case for any signs of damage or wear before each use.
- Keep the case out of reach of children and vulnerable individuals when not in use.
- Follow all instructions provided in this guide to avoid potential hazards.
- In case of an emergency, familiarize yourself with the rapid deployment of the tourniquet.

Specific Safety Precautions for Use

- **Avoid Excessive Force:** Do not apply excessive force when pulling the tourniquet from the case, as this may damage the case or the tourniquet.
- **Moisture Management:** Ensure that the case and tourniquet are kept dry to maintain functionality. If wet, dry thoroughly before use.
- **Secure Mounting:** Ensure that the case is securely mounted to your duty belt, vest, or holster to prevent accidental detachment.
- **Regular Maintenance:** Clean the case regularly to remove dirt and debris that may affect its performance.
- **Training:** Only individuals trained in the use of a tourniquet should attempt to deploy it from this case.

Instructions for Installation and Usage

1. Mounting the Case:

- Choose a suitable location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- Use the appropriate mounting hardware (sold separately for some options) to secure the case in place.
- Ensure that the mounting is tight and stable to prevent movement during use.

2. Loading the Tourniquet:

- Insert the CAT® Tourniquet (Gen 7 or older) into the case, ensuring it is seated securely.
- The tourniquet should be positioned for easy access, with the exposed portion facing outward.

3. Deploying the Tourniquet:

- In an emergency, grasp the exposed portion of the tourniquet firmly.
- Pull swiftly to release the tourniquet from the case.
- Follow the manufacturer's instructions for applying the tourniquet effectively.

4. Adjusting the Height:

- If necessary, adjust the mounting position of the case to ensure it is at a comfortable height for

quick access.

- The case can be mounted in vertical, horizontal, canted, or inverted positions as needed.

Disposal Instructions

- Dispose of the Rigid TQ Case in accordance with local regulations for plastic waste.
- If the case is damaged or no longer functional, ensure it is disposed of properly to prevent environmental harm.

Contact Information for Further Support

For any safety inquiries or concerns regarding the Rigid TQ Case, please reach out to the appropriate support channels as indicated by the retailer or manufacturer.

Thank you for your attention to these safety guidelines. Your safety and the safety of others depend on the proper use of this product.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk