

# RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



## Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, COYOTE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043726
- Mfr. No.: E10-3011M-CYT
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668274

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions](#)
- [About Us](#)

# RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions

## Introduction

Thank you for choosing the RIGID TQ Case for SOF®TT/SOF®TTW. This safety instruction guide is designed to help you use this product safely and effectively. The RIGID TQ Case is a hardsided tourniquet case that ensures rapid access to your tourniquet when every second counts. Please read this guide carefully to understand how to use and care for your product.

## General Safety Guidelines

- Ensure that the RIGID TQ Case is used only for its intended purpose, which is to hold SOF®TT or SOF®TTWide tourniquets.
- Regularly inspect the case for any signs of wear or damage, which may affect its performance.
- Keep the case away from extreme temperatures, moisture, and direct sunlight to maintain its integrity.
- Store the RIGID TQ Case in a secure location to prevent unauthorized access, especially around children.
- If you notice any defects or issues with the product, discontinue use immediately and seek further assistance.

## Specific Safety Precautions for Use

- Always check that the tourniquet is securely held in place within the RIGID TQ Case before use.
- Use the case as instructed to ensure rapid deployment of the tourniquet in an emergency.
- Avoid using excessive force when removing the tourniquet from the case, as this may damage the case or the tourniquet.
- Ensure that the mounting solutions used (e.g., MOLLE/PALS, drop leg holster) are securely attached and compatible with your equipment.
- If using the case while engaged in strenuous activities, ensure that it is properly secured to prevent accidental detachment.

## Instructions for Installation and Usage

### Installation

#### 1. Mounting Options:

- The RIGID TQ Case can be mounted on:
  - Police Duty Belt using the BladeTech TEKLOK Belt Attachment (fits up to a 2.25" belt).
  - MOLLE/PALS compatible vest.
  - Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately).

#### 2. Securing the Case:

- Use the included Vibratite to secure the mounting hardware once the desired configuration is obtained.
- Follow the instructions provided with the Vibratite for proper application.

### Usage

#### 1. Accessing the Tourniquet:

- To deploy the tourniquet, grasp the exposed portion firmly.
- With a swift pull, the tourniquet will release from the case.
- Ensure that your hands are in a position to handle the tourniquet effectively, even if they are wet or bloody.

## **2. Maintenance:**

- Clean the RIGID TQ Case with a damp cloth as needed. Do not use harsh chemicals that may damage the material.
- Regularly check the drainage hole at the bottom of the case to ensure it is clear and functioning.

## **Disposal Instructions**

- Dispose of the RIGID TQ Case according to local regulations for plastic waste.
- If the case is damaged beyond repair, consider recycling the materials if possible.

## **Contact Information for Further Support**

For any questions or concerns regarding the RIGID TQ Case, please refer to the manufacturer's contact details provided on the packaging or the official website.

## **Conclusion**

The RIGID TQ Case for SOF®TT/SOF®TTW is designed to provide secure storage and rapid access to your tourniquet. By following these safety instructions, you can ensure that you are prepared for emergencies while using this product effectively. Always prioritize safety and stay informed about best practices for using your RIGID TQ Case.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)