

RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, ORANGE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, ORANGE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043724
- Mfr. No.: E10-3011B-ORG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668380

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions](#)
- [About Us](#)

RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions

Introduction

Thank you for choosing the RIGID TQ Case for SOF®TT/SOF®TTW. This product is designed to provide secure storage and rapid access to your tourniquet in emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the RIGID TQ Case is used only for its intended purpose, which is to hold SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of damage or wear. If any issues are found, discontinue use immediately.
- Keep the case and its contents out of reach of children and vulnerable individuals to prevent accidental use.
- Always ensure that the tourniquet is properly secured within the case to avoid accidental release.
- Follow all local regulations and guidelines regarding the use of medical devices.

Specific Safety Precautions for Use

- Use the RIGID TQ Case only with compatible tourniquets (SOF®TT or SOF®TTW) to ensure proper fit and function.
- Avoid exposing the case to extreme temperatures, chemicals, or harsh environments that may compromise its integrity.
- When deploying the tourniquet, ensure that you have a clear line of sight and space to operate effectively.
- If the case becomes wet or contaminated, clean it according to the manufacturer's instructions to maintain hygiene and functionality.
- Always practice deploying the tourniquet in a safe environment to familiarize yourself with the process before an emergency arises.

Instructions for Installation and Usage

Installation

- Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- Use the provided mounting options:
 - For Police Duty Belt: Attach using the BladeTech TEKLOK Belt Attachment, which fits belts up to 2.25 inches.
 - For MOLLE/PALS Attachment: Use the Tactical Tailor Malice Clip.
- If mounting to a Safariland® Drop Leg Shroud, ensure that the accessory position is used, and obtain any necessary mounting hardware.
- Utilize the included Vibratite to secure the hardware in place once the desired configuration is achieved.

Usage

- To deploy the tourniquet:
 - Grasp the exposed portion of the tourniquet firmly.
 - With a swift pull, release the tourniquet from the case.
- Ensure that you are familiar with the tourniquet application process and practice as needed.

Disposal Instructions

- Dispose of the RIGID TQ Case in accordance with local regulations regarding plastic and medical waste.
- If the case is damaged beyond repair, ensure it is disposed of safely to prevent any potential harm.

Contact Information for Further Support

For support regarding the RIGID TQ Case, please refer to the manufacturer's website or contact the retailer from whom you purchased the product for assistance.

By following these safety instructions, you can ensure the effective use of the RIGID TQ Case for SOF®TT/SOF®TTW and maintain a safe environment for yourself and others. Always stay informed about product recalls and safety alerts through the EU's Safety Gate platform. Thank you for prioritizing safety.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk