

# RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, RANGER GREEN

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



## Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, RANGER GREEN](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043723
- Mfr. No.: E10-3011B-RGR
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668298

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions](#)
- [About Us](#)

# RIGID TQ Case for SOF®TT/SOF®TTW Safety Instructions

## Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This product is designed to securely hold your tourniquet, providing quick access when it matters most. To ensure safe and effective use of this product, please read and follow all safety instructions outlined in this guide.

## General Safety Guidelines

- **Product Safety:** This product is designed for safe use. Ensure that you understand how to use it correctly to avoid any potential hazards.
- **Enhanced Recalls:** In the event of a safety issue, recall notices will be provided. Always check for updates regarding product safety.
- **Online Shopping:** If purchased online, ensure the seller complies with safety regulations similar to those in physical stores.
- **Special Consumer Focus:** Extra precautions should be taken when using this product around vulnerable individuals, including children.
- **EU Contact Point:** For any safety inquiries, please refer to the designated EU contact point.
- **Rapid Alerts:** Stay informed about product safety by checking the EU's Safety Gate platform for updates on unsafe products.

## Specific Safety Precautions for Use

- **Hazard Identification:**
  - Ensure that the tourniquet is securely held in the case to prevent accidental release.
  - Be cautious of sharp edges or components during installation and use.
- **Avoiding Hazards:**
  - Do not use the case if it is damaged or shows signs of wear.
  - Ensure that the mounting solution is securely attached to your belt or vest to prevent it from falling.
  - Avoid exposing the product to extreme temperatures or harsh chemicals, as this may affect its integrity.
- **Age Specific Warnings:**
  - This product is not intended for use by children. Always keep it out of reach of minors.

## Instructions for Installation and Usage

- **Installation:**
  1. Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
  2. If using a BladeTech TEKLOK Belt Attachment, ensure it fits securely on belts up to 2.25 inches.
  3. For MOLLE/PALS attachment, use the Tactical Tailor Malice Clip to secure the case.
  4. If mounting to a Safariland® Drop Leg Shroud, ensure the mounting hardware is obtained separately.
  5. Use the included Vibratite to secure all hardware in place once the desired configuration is obtained.
- **Usage:**
  1. To deploy the tourniquet, grasp the exposed portion firmly.

2. With a swift pull, the tourniquet will release from the case.
3. Practice using the case to ensure you can deploy the tourniquet effectively under pressure.
4. Regularly check the case and tourniquet for any signs of wear or damage.

## **Disposal Instructions**

- Dispose of this product in accordance with local regulations.
- If the product is no longer needed, ensure that it is disposed of in a way that prevents accidental use by others, especially children.

## **Contact Information for Further Support**

For further assistance or inquiries regarding the RIGID TQ Case for the SOF®TT/SOF®TTW, please refer to the manufacturer's contact information provided with your product packaging or on their official website.

Thank you for your attention to these safety instructions. Your safety and the safety of others depend on the proper use of this product.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)